


TECHNICAL DATA SHEET

ANTIX HIP PROTECTOR



- The protector is used to protect the reproductive organs.
- Level protection: front Pb mm 0.50.
- On request: embroidery - name department.
- Textile and edging of your choice for easy identification in the department.
- CE, ISO9001 quality system and compliance to EN 61331-3.

CODE	DESCRIPTION
STANDARD MODEL WEIGHT: 3,8 kg mq Radioprotective classic material lead-based and a film of barium in surface	
O-GB-70	Kids – Dim. cm 20 x 25
O-GB-71	Kids – Dim. cm 30 x 30
O-GB-72	Adults – Dim. cm 40 x 40
O-GB-73	Adults – Dim. cm 50 x 50