

**SCHERMATURE ANTIX - ACCESSORI PER RADIOLOGIA****CE CONFORMITY DECLARATION**

The undersigned "CABLAS S.r.l."  
office at Via Michelangelo Buonarroti 52/B  
20093 Cologno Monzese (MI) - ITALY

DECLARES

Under own responsibility that the product

**MOBILE SCREENS WITH ANTIX LEAD GLASS**

Codes

S-SM-20	S-SM-21	S-SM-22	S-SM-23
---------	---------	---------	---------

S-SM-24	S-SM-25	S-SM-26	S-SM-27
---------	---------	---------	---------

This declaration is conforms to the norm "GOOD TECHNICAL"

As required by Art. 72 D. Decree 230 of 1995, which transposed the Directive EURATOM 80/836, 84/466, 89/618, 90/641, 92/3 in the field of ionizing radiation.

They are produced in our specialized workshop and emplaced a "workmanlike".


Lead is used first melting purposes only 99.9% in conformity with our UNI 3165,  
which corresponds to 2 mm Pb

The crystal antix complies with the test report no. 49046 made by the manufacturer at specialized laboratories. The technical features are related to the IEC 1331-1 (1994-10)

**TUBE VOLTAGE RX: 150KV - GLASS THICKNESS: MM7 / 8.5**

Cologno Monzese, year 2015

**CABLAS S.R.L.**  
Via Michelangelo, 52/B  
20093 COLOGNO MONZESE (MI)  
S.r.l.: Via Sormani 34 - COLOGNO M.SE  
Tel. R.A. 02 2538088 - Fax 02 27300644  
Cod. Fisc. e P. IVA 11903230156



---

**SCHERMATURE ANTIX - ACCESSORI PER RADIOLOGIA**

---