

## TECHNICAL DATA - SHEET

### ANTIX APRON COAT




*The product  
don't require the  
GMDN code  
as  
**INDIVIDUAL  
ANTIX  
PROTECTION  
DEVICE***

- Model double-breasted, provides protection to 360 degrees.
- The two front panels double the protection.
- Allows you to sit and bend easily.
- Adjustable belt in height.
- Textile and edging of your choice for easy identification in the department.
- Sizes shaped in accordance with EN 340
- CE, ISO9001 quality system and compliance to EN 61331-3
- Pocket and padded shoulder straps always included

## TECHNICAL DATA - SHEET

CODE	SIZE	FRONTAL PROTECTION	BACK PROTECTION	LENGHT
<b>STANDARD MODEL WEIGHT: 3,8Kg Mq</b>				
Radioprotective classic material lead-based and a film of barium in surface				
<b>NC25100</b>	S / M / L	PB MM 0,25	PB MM 0,25	CM 100
<b>NC25110</b>	S / M / L	PB MM 0,25	PB MM 0,25	CM 110
<b>NC25120</b>	S / M / L	PB MM 0,25	PB MM 0,25	CM 120
<b>NC50100</b>	S / M / L	PB MM 0,50	PB MM 0,25	CM 100
<b>NC50110</b>	S / M / L	PB MM 0,50	PB MM 0,25	CM 110
<b>NC50120</b>	S / M / L	PB MM 0,50	PB MM 0,25	CM 120
<b>MODEL LIGHTENED WEIGHT: 3,2 Kg Mq</b>				
Radioprotective material with the addition of metals which replace partially lead, reducing weight				
<b>NC25100/A</b>	S / M / L	PB MM 0,25	PB MM 0,25	CM 100
<b>NC25110/A</b>	S / M / L	PB MM 0,25	PB MM 0,25	CM 110
<b>NC25120/A</b>	S / M / L	PB MM 0,25	PB MM 0,25	CM 120
<b>NC50100/A</b>	S / M / L	PB MM 0,50	PB MM 0,25	CM 100
<b>NC50110/A</b>	S / M / L	PB MM 0,50	PB MM 0,25	CM 110
<b>NC50120/A</b>	S / M / L	PB MM 0,50	PB MM 0,25	CM 120
<b>MODEL "NO LEAD" WEIGHT: 2,7Kg Mq</b>				
High-tech material that replaces the lead with metals which are balanced from the point of view of atomic				
<b>NC25100/NL</b>	S / M / L	PB MM 0,25	PB MM 0,25	CM 100
<b>NC25110/NL</b>	S / M / L	PB MM 0,25	PB MM 0,25	CM 110
<b>NC25120/NL</b>	S / M / L	PB MM 0,25	PB MM 0,25	CM 120
<b>NC50100/NL</b>	S / M / L	PB MM 0,50	PB MM 0,25	CM 100
<b>NC50110/NL</b>	S / M / L	PB MM 0,50	PB MM 0,25	CM 110
<b>NC50120/NL</b>	S / M / L	PB MM 0,50	PB MM 0,25	CM 120



<b>CORRESPONDENCE SIZE "WOMEN"</b>			
Based on the norm EN340.2004			
	<b>S</b>	<b>M</b>	<b>L</b>
Circumf. breast cm	94	98	102
Waist circumf. cm	74	78	82
Hip circumf.	100	104	108

<b>CORRESPONDENCE SIZE "MAN"</b>			
Based on the norm EN340.2004			
	<b>S</b>	<b>M</b>	<b>L</b>
Thorax circumf. cm	92	96	100
Waist circumf. cm	80	84	88

