



Via Michelangelo Buonarroti, 52/B - 20093 Cologno M.se (Mi)

## TECHNICAL DATA-SHEET

### INTENSIFYNG SCREEN RARE EARTH BLUE / GREEN GREEN GRADUAL



The intensifying screens are compatible with all conventional x-ray films and are destined to improve the quality of overall rx image.

Our range of screens is suitable for all type of x-ray cassettes.

Every screen provides high contrasts, details perception, even with low rx radiation dosage.

## DIMENSION

### INTENSIFYING SCREEN RARE EARTH BLUE - SPEED 400

CODE	SIZE
H-AS-01	13 X 18
H-AS-02	18 X 24
H-AS-03	24 X 30
H-AS-04	18 X 43
H-AS-05	30 X 40
H-AS-06	35 X 35
H-AS-07	35 X 43
H-AS-08	15 X 30
H-AS-09	15 X 40
H-AS-10	20 X 40
H-AS-12	30 X 90
H-AS-13	30 X 120

### INTENSIFYING SCREEN RARE EARTH GREEN GRADUAL

CODE	SIZE
H-GR-09	15 X 40
H-GR-10	20 X 40
H-GR-12	30 X 90
H-GR-13	30 X 120
H-GR-14	40 X 90
H-GR-15	40 X 120

**INTENSIFYING SCREEN RARE EARTH GREEN - SPEED 400**

<b>CODE</b>	<b>SIZE</b>
H-AT-20	13 X 18
H-AT-21	18 X 24
H-AT-22	35 X 43
H-AT-23	18 X 43
H-AT-24	30 X 40
H-AT-25	35 X 35
H-AT-26	35 X 43
H-AT-27	15 X 30
H-AT-28	15 X 40
H-AT-29	20 X 40
H-AT-31	30 X 90
H-AT-32	30 X 120