

TECHNICAL DATA - SHEET

NON-MAGNETIC CHAIR



- Tubular aluminium frame;
- Seat and back upholstered with soft washable synthetic leather, folding footrest, capacity 135 kg;
- Sliding on No. 4 swivel wheels stainless steel 18/10 mm diameter 100 of which No. 2 with brakes;
- Handling by No. 2 comfortable non-slip rubber grips.

CODE	DESCRIPTION
R-MN-10	Outside Dim. cm 55 x 70 x 90 h Max weight carried Kg 135 Weight Kg 16.