

TECHNICAL DATA-SHEET

MOBILE SCREENS WITH ANTIX LEAD GLASS



- Provided with x ray panel from in the lower part with a plastic laminate coating on both sides.
- Antix visual glass in the upper part over the entire width.
- Screening Pb 2 mm.
- Outside aluminium structure.
- Rubbered wheels with brakes.

CODE	SHIELDING	OUTSIDE DIMENSION	VISUAL	WEIGHT
S-SM-20	Pb mm 2	CM 82 X 198H	CM 80 X 60H	KG 100
S-SM-21	Pb mm 2	CM 82 X 198H	CM 80 X 80H	KG 98
S-SM-22	Pb mm 2	CM 102 X 198H	CM 100 X 60H	KG 112
S-SM-23	Pb mm 2	CM 102 X 198H	CM 100 X 80H	KG 120
S-SM-24	Pb mm 2	CM 122 X 198H	CM 120 X 60H	KG 130
S-SM-25	Pb mm 2	CM 122 X 198H	CM 120 X 80H	KG 136
S-SM-26	Pb mm 2	CM 142 X 198H	CM 140 X 60H	KG 146
S-SM-27	Pb mm 2	CM 142 X 198H	CM 140 X 80H	KG 150