

## TECHNICAL DATA - SHEET VISORE SPOT

- Suitable for a reading of x-rays overlapped.
- Structure in painted steel sheet in focus.
- Alogen high intensity light lamp 300 W, long life.
- Suitable for exposed films reading.



CODE	DESCRIPTION
	<b>TABLE SPOT</b>
<b>M-SP-01</b>	<ul style="list-style-type: none"><li>• With base</li><li>• push-button switch</li></ul>
	<b>WALL SPOT</b>
<b>M-SP-02</b>	<ul style="list-style-type: none"><li>• Equipped with support and pole cm43h</li><li>• push-button switch</li></ul>
	<b>ILLUMINATOR SPOT</b>
<b>M-SP-03</b>	<ul style="list-style-type: none"><li>• Equipped with support and pole for outside application</li><li>• push-button switch</li></ul>
	<b>TABLE SPOT</b>
<b>M-SP-05</b>	<ul style="list-style-type: none"><li>• With base</li><li>• switch with pedal.</li></ul>

### ACCESSORIES

<b>M-DI-01</b>	IRIS DIAPHRAGM
	<b>LENS ON FLEXIBLE ARM</b>
<b>M-LF-01</b>	<ul style="list-style-type: none"><li>• Diameter 120 mm, complete with attack spot.</li></ul>