


TECHNICAL DATA - SHEET

X – RAY HANGING PROTECTIONS



CODE	DESCRIPTION
V-SP-11	Stand to the ceiling with glass 50 x 40
V-SP-23	Surgical lamp 35.000 LUX RA 93/1 mt with fixed focus Complete with movable arm electrified.
V-SP-24	Surgical lamp 50.000 LUX RA 96/1mt With adjustable light 14-25. Complete with movable arm electrified.
V-SP-32	Castle Extension for ceiling anchors. Suitable for heights "floor to ceiling" greater than 3 m
V-SP-33	Sliding system ceiling Maximum length: 250 cm

- For eyes and thyroid protection to be used in Angiography, urology, dpts etc.
- Provided with a ceiling support by a mobile arm (length cm 80) articulated and balanced (it can turn 325 deg.) and by a further arm (length cm 80) moving up and down 50 deg.(it can turn 360 deg)
- This item has at the lowest extremity a frame with glass X-ray Pb mm 0.50 2.0 50 x 40 cm , and one little apron X- ray Pb mm 0.50 (it can turn 360 degree)
- Prepared for heights floor / ceiling Mt from 2.70 mt to 3.00, to be specified at time of order