

TECHNICAL DATA - SHEET

X- RAY PROTECTION VASCULAR RADIOLOGY

V-SP-41



- Consists of a metal structure that supports the top No. 1 mm Pb acrylic crystal antix 0:50 shaped 116x50h cm;
- In the lower part No. 2 towels antix 48 x 80 h - Pb 0.50 mm which one extendable horizontally;
- Both are adjustable in height by means of the balancing system 146 x 177 cm;;
- Suitable for use in radiological vascular surgery..

CODE	DESCRIPTION
V-SP-41	Screen protection for vascular surgery with shaped crystal.