

TECHNICAL DATA - SHEET


ANTIX APRON WITH SCAPULAR PROTECTION SINGLE SIZE



- Complete frontal protection.
- Padded shoulders for reduced shoulder stress and equitable distribution of weight.
- Belt snap closure.
- Easy to wear and remove.
- Pocket included.
- Allows you to sit and bend easily.

TECHNICAL DATA-SHEET

CODE	FRONTAL PROTECTION	BACK PROTECTION	WIDTH	LENGHT
LIGHT MODEL				
KGB10/100	PB MM 0,25	PB MM 0,25	CM 60	CM 100
KGB10/110	PB MM 0,25	PB MM 0,25	CM 60	CM 110
KGB11	PB MM 0,25	PB MM 0,25	CM 60	CM 120
KGB07/100	PB MM 0,35	PB MM 0,25	CM 60	CM 100
KGB07/110	PB MM 0,35	PB MM 0,25	CM 60	CM 110
KGB08	PB MM 0,35	PB MM 0,25	CM 60	CM 120
KGB12/100	PB MM 0,50	PB MM 0,25	CM 60	CM 100
KGB12/110	PB MM 0,50	PB MM 0,25	CM 60	CM 110
KGB13	PB MM 0,50	PB MM 0,25	CM 60	CM 120



CORRESPONDENCE SIZE			
"WOMEN- MAN"			
<small>Based on the norm EN340.2004</small>			
	SINGLE SIZE		
Lenght	100	110	120
Width	60		

