


TECHNICAL DATA - SHEET

ANTIX COLLARS WITH OR WITHOUT BIB



- This model protects the front of the neck (thyroid) and fastens with Velcro in the back side;
- Coated with elegant teflon nylon, water, oil and dirt repellent;
- Any impurities stay on the surface and can be easily cleaned with a soft towel;
- Is appropriate to use it with the antix apron;
- CE, ISO9001 quality system and compliance to EN 61331-3

CODICE	NECK	PROTECTION	BIB
STANDARD MODEL WEIGHT: 3,8Kg Mq			
Radioprotective classic material lead-based and a film of barium in surface			
O-GD-35	CM 29/40	PB MM 0,25	NO
O-GD-36	CM 29/40	PB MM 0,50	NO
O-GD-37	CM 29/40	PB MM 0,25	YES
O-GD-38	CM 29/40	PB MM 0,35	YES
O-GD-40	CM 29/40	PB MM 0,50	YES
CODE	NECK	PROTECTION	HARNESS
MODEL "NO LEAD" WEIGHT: 2,7Kg Mq			
High-tech material that replaces the lead with metals which are balanced from the point of view of atomic			
O-GD-50	CM 29/40	PB MM 0,25	YES
O-GD-51	CM 29/40	PB MM 0,35	YES
O-GD-52	CM 29/40	PB MM 0,50	YES