



Via Michelangelo Buonarroti, 52/B - 20093 Cologno M.se (Mi)

## TECHNICAL DATA SHEET

### **RADIOGRAPHIC CASSETTES FOR MAMMOGRAPHY WITH WINDOW**



The material with which the cassette has been built, molded polycarbonate, allows a quick and easy loading of films, both manual and automatic.

The single mechanical lock makes the box sturdier and more durable.

The box allows you to place the film in an extremely close to the chest wall of the patient, ensuring excellent quality.

**DIMENSION**  
(COMPLETE SCREENS OF REINFORCED MONO - EMULSION)

<b>CODE</b>	<b>SIZE</b>
G-MM-03	18 X 24
G-MM-04	24 X 30