


CABLAS s.r.l.

Via Michelangelo Buonarroti, 52/B - 20093 Cologno M.se (Mi)

TECHNICAL DATA-SHEET

RADIOGRAPHIC CASSETTES



The x-ray cassettes are produced in alloy aluminium and special material, the edges are in plastic and aluminium.

The internal covering material is realized with a soft synthetic fabric that allows the intensifying screen to stick perfectly to the x-ray film, assuring that the image on the film is clear and rich.

DIMENSIONS OF RADIOGRAPHIC CASSETTES

| CODE | SIZE |
|----------------|-------------|
| G-SP-50 | 13x18 |
| G-SP-51 | 18x24 |
| G-SP-52 | 24x30 |
| G-SP-54 | 30x40 |
| G-SP-55 | 35x35 |
| G-SP-56 | 35x43 |
| G-SP-57 | 15x30 |
| G-SP-58 | 15x40 |
| G-SP-59 | 20x40 |
| G-SP-60 | 30x80 |
| G-SP-61 | 30x90 |
| G-SP-62 | 30x120 |
| G-SP-63 | 40x120 |
| G-SP-67 | 40x90 |
| G-SP-69 | 18x43 |