

TECHNICAL DATA - SHEET

X-RAY SHEETS



- Suitable for shielding of windows and doors;
- Made of extraflexible leaded polyvnyl sheets ;
- Nylon covered with perimeter edges

CODE	DESCRIPTION
	STANDARD MODEL WEIGHT: 3,8Kg Mq Radioprotective classic material lead-based and a film of barium in surface
O-TA-01	Antix curtain Pb mm 0,25
O-TA-02	Antix curtain Pb mm 0,50
O-TA-50	Curtain track