


## TECHNICAL DATA-SHEET

### ANTIDIFFUSION X-RAY GRIDS 200 L/P FOR MAMMOGRAPHY FOCUSED CM 50/60



This kind of x-ray grids support professional doctors to make accurate diagnosis improving the contrast of the image, absorbing scattered x-ray and maintaining an optimal transmission of primary radiation:

#### **A**

Contrast: The antidiffusion x-ray grids are produced using products of high purity

#### **B**

Ingredients: The x-ray grids use lead and aluminium of high quality which allow to obtain an image with such a high contrast, that you can see the difference with your naked eyes.

**DIMENSION**

**ANTIDIFFUSION X-RAY GRIDS 200 L/ P  
FOR MAMMOGRAPHY FOCUSED CM 50/60**

**RATIO AVAILABLE: 5:1**

<b>CODE</b>	<b>SIZE</b>	<b>WEIGHT Kg</b>
<b>I-GR-420</b>	18 X 24	0.486
<b>I-GR-421</b>	24 X 30	0.724