

TECHNICAL DATA - SHEET

ANTIX SKIRT AND VEST LONG PROCEDURES



*The product
don't require the
GMDN code
as
**INDIVIDUAL
ANTIX
PROTECTION
DEVICE***



- Model double-breasted, provides protection to 360 degrees.
- Allows you to sit and bend easily.
- Adjustable belt in height.
- Pocket and padded shoulder straps always included.
- Skirt and full overlay front.
- Completely enveloping
- The front panels double the protection.
- Non-slip rear interior rear.
- Long slots support.
- Double closure: VELCRO + BELT.
- Sizes shaped man / woman in compliance with EN 340, CE Mark, ISO9001 quality system and compliance to EN 61331-3.

TECHNICAL DATA - SHEET

CODE	SIZE	FRONTAL PROTECTION	BACK PROTECION
STANDARD MODEL WEIGHT: 3,8 kg mq			
Radioprotective classic material lead-based and a film of barium in surface			
N2PS25	S / M / L	PB MM 0,25	PB MM 0,25
N2PS35	S / M / L	PB MM 0,35	PB MM 0,25
N2PS50	S / M / L	PB MM 0,50	PB MM 0,25
MODEL LIGHTENED WEIGHT: 3,2 Kg Mq			
Radioprotective material with the addition of metals which replace partially lead, reducing weight			
N2PS25/A	S / M / L	PB MM 0,25	PB MM 0,25
N2PS35/A	S / M / L	PB MM 0,35	PB MM 0,25
N2PS50/A	S / M / L	PB MM 0,50	PB MM 0,25
MODEL "NO LEAD" WEIGHT: 2,7Kg Mq			
High-tech material that replaces the lead with metals which are balanced from the point of view of atomic			
N2PS25/NL	S / M / L	PB MM 0,25	PB MM 0,25
N2PS35/NL	S / M / L	PB MM 0,35	PB MM 0,25
N2PS50/NL	S / M / L	PB MM 0,50	PB MM 0,25

CORRESPONDENCE SIZE "WOMEN"			
Based on the norm EN340.2004			
	S	M	L
Circumf. breast cm	94	98	102
Waist circumf. cm	74	78	82
Hip circumf.	100	104	108

CORRESPONDENCE SIZE "MAN"			
Based on the norm EN340.2004			
	S	M	L
Thorax circumf.cm	92	96	100
Waist circumf. cm	80	84	88

