

TECHNICAL DATA - SHEET

ANTIX BLANKET



The antix blanket is used for the screening of patients.

- Made of flexible antix rubber.
- Coated water-repellent nylon blue color.
- Perimeter edging tear.
- CE, ISO9001 quality system and compliance to EN 61331-3.

CODE	DESCRIPTION
O-TA-60	Antix blanket Pb mm 0,50 – Dim. cm 100 x 80
O-TA-61	Antix blanket Pb mm 0,50 – Dim. cm 120 x 80
O-TA-62	Antix blanket Pb mm 0,50 – Dim. cm 150 x 80
O-TA-63	Antix blanket Pb mm 0,50 – Dim. cm 100 x 100
O-TA-64	Antix blanket Pb mm 0,50 – Dim. cm 120 x 100
O-TA-65	Antix blanket Pb mm 0,50 – Dim. cm 150 x 100