

TECHNICAL DATA - SHEET

SCROTAL AND OVARIES PROTECTIONS



Suitable for the protection of the genital organs.

CODE	DESCRIPTION
O-GI-05	SCROTAL PROTECTIONS PB MM 2 Complete series n° 4 sizes
O-GI-22	OVARIES PROTECTIONS PB MM 1 Complete series n° 5 sizes: • PLATE 1: cm16 x 15.5 • PLATE 2: cm 15.5 x 15 • PLATE 3: cm 13.5 x 14.5 • PLATE 4: cm 12.5X 13.5 • PLATE 5: cm 11 x 10.5