

## TECHNICAL DATA - SHEET


### SURGERY ANTIX APRON



*The product  
don't require the  
GMDN code  
as  
**INDIVIDUAL  
ANTIX  
PROTECTION  
DEVICE***

- It is removable from under a sterile gown.
- Detaching the bands to tear the gown slips down.
- Kept here almost exclusively in surgery.
- The model was shaped to allow large movements of the operator.
- Three levels protection: front Pb mm 0.25 / 0.35 / 0.50mm – back Pb mm 0.25.
- Available also without back protection.
- Textile and edging of your choice for easy identification in the department.
- Sizes shaped in accordance with EN 340.
- CE, ISO9001 quality system and compliance to EN 61331-3.
- On request: embroidery-name department / elastic belt.
- Pocket and padded shoulder straps always included.

CODE	SIZE	FRONTAL PROTECTION	BACK PROTECTION	LENGHT
<b>STANDARD MODEL WEIGHT: 3,8Kg Mq</b>				
Radioprotective classic material lead-based and a film of barium in surface				
NCH25100	S / M / L	PB MM 0,25	PB MM 0,25	CM 100
NCH25110	S / M / L	PB MM 0,25	PB MM 0,25	CM 110
NCH25120	S / M / L	PB MM 0,25	PB MM 0,25	CM 120
NCH35100	S / M / L	PB MM 0,35	PB MM 0,25	CM 100
NCH35110	S / M / L	PB MM 0,35	PB MM 0,25	CM 110
NCH35120	S / M / L	PB MM 0,35	PB MM 0,25	CM 120
NCH50100	S / M / L	PB MM 0,50	PB MM 0,25	CM 100
NCH50110	S / M / L	PB MM 0,50	PB MM 0,25	CM 110
NCH50120	S / M / L	PB MM 0,50	PB MM 0,25	CM 120
<b>MODEL LIGHTENED WEIGHT: 3,2 Kg Mq</b>				
Radioprotective material with the addition of metals which replace partially lead, reducing weight				
NCH25100/A	S / M / L	PB MM 0,25	PB MM 0,25	CM 100
NCH25110/A	S / M / L	PB MM 0,25	PB MM 0,25	CM 110
NCH25120/A	S / M / L	PB MM 0,25	PB MM 0,25	CM 120
NCH35100/A	S / M / L	PB MM 0,35	PB MM 0,25	CM 100
NCH35110/A	S / M / L	PB MM 0,35	PB MM 0,25	CM 110
NCH35120/A	S / M / L	PB MM 0,35	PB MM 0,25	CM 120
NCH50100/A	S / M / L	PB MM 0,50	PB MM 0,25	CM 100
NCH50110/A	S / M / L	PB MM 0,50	PB MM 0,25	CM 110
NCH50120/A	S / M / L	PB MM 0,50	PB MM 0,25	CM 120
<b>MODEL "NO LEAD" WEIGHT: 2,7 Kg Mq</b>				
High-tech material that replaces the lead with metals which are balanced from the point of view of atomic				
NCH25100/NL	S / M / L	PB MM 0,25	PB MM 0,25	CM 100
NCH25110/NL	S / M / L	PB MM 0,25	PB MM 0,25	CM 110
NCH25120/NL	S / M / L	PB MM 0,25	PB MM 0,25	CM 120
NCH35100/NL	S / M / L	PB MM 0,35	PB MM 0,25	CM 100
NCH35110/NL	S / M / L	PB MM 0,35	PB MM 0,25	CM 110
NCH35120/NL	S / M / L	PB MM 0,35	PB MM 0,25	CM 120
NCH50100/NL	S / M / L	PB MM 0,50	PB MM 0,25	CM 100
NCH50110/NL	S / M / L	PB MM 0,50	PB MM 0,25	CM 110
NCH50120/NL	S / M / L	PB MM 0,50	PB MM 0,25	CM 120



<b>CORRESPONDENCE SIZE "WOMEN"</b>			
Based on the norm EN340.2004			
	S	M	L
Circumf. breast cm	94	98	102
Waist circumf. cm	74	78	82
Hip circumf.	100	104	108



<b>CORRESPONDENCE SIZE "MAN"</b>			
Based on the norm EN340.2004			
	S	M	L
Thorax circumf. cm	92	96	100
Waist circumf. cm	80	84	88