

## TECHNICAL DATA-SHEET

### SLIM X -RAY FILM VIEWE



	<b>COD. LED1 SECTION 1</b>	<b>COD. LED2 SECTION 2</b>	<b>COD. LED3 SECTION 3</b>	<b>COD. LED4 SECTION 4</b>
<b>Net Weight</b>	6 kg	9 kg	16 kg	21.5 kg
<b>Gross Weight</b>	8 kg	11 kg	18.5 kg	25.5 kg
<b>Outside dimensions</b>	504x502x25 mm	884x502x25 mm	1265x502x25 mm	1646x502x25 mm
<b>Dimension visual field</b>	362x440 mm	742x440 mm	1123x440 mm	1504x440 mm
<b>Thickness</b>	25mm	25mm	25mm	25mm
<b>N. of lamps</b>	4 lamps	8 lamps	12 lamps	16 lamps
<b>Lamp life</b>	100,000hours			
<b>Models lamps</b>	LED			
<b>Color Temperature</b>	≥ 6500k			
<b>Luminosity</b>	Adjustable 0-2000cd/m <sup>2</sup> (6500Lux)			
<b>Uniformity</b>	≥ 98%			
<b>Frequency</b>	50KHz			
<b>Voltage</b>	90V-260V, 50Hz±1Hz			
<b>Suitability films</b>	Normal films, digital high-intensity, mammography films.			
<b>Might</b>	15w	30w	45w	60w

