






## TRCHNICAL DATA - SHEET

### TANKS OF DEVELOPMENT MANUAL

<b>A-VS-15</b> KIT DEVELOPMENT FOR FRAMES	<b>A-VS-20</b> ALFA 4	<b>A- VS-21</b> ALFA 4 THERMOREGULATED	<b>A-VS-22</b> TIPO A	<b>A-VS-24</b> TIPO A THERMOREGULATED
				
<p>Manual development of films rx up to the FORMAT CM 35 X 43</p> <p><b>INCLUDED:</b></p> <ul style="list-style-type: none"> <li>-2 20 lt containers with drain cock</li> <li>-1 30 lt containers with drain cock.</li> <li>-1 Tank for containing liquid drain containers</li> <li>-5 M rubber tube flexible.</li> </ul> <p><b>DIMENSION :</b> Cm 50x60x62</p> <p><b>WEIGHT:</b> 22 Kg</p>	<p>Manual development of films rx up to the FORMAT CM 35 X 43</p> <p><b>INCLUDED:</b></p> <ul style="list-style-type: none"> <li>-2 fixed compartments by 5 liters for developing and fixing</li> <li>-1 Fixed compartment by 28 lt for intermediate and final washing</li> <li>- Automatic loading and overflow</li> </ul> <p><b>DIMENSION :</b> Cm 27 x 48 x 96</p> <p><b>WEIGHT:</b>12 Kg</p>	<p>Manual development of films rx up to the FORMAT CM 35 X 43</p> <p><b>INCLUDED:</b></p> <ul style="list-style-type: none"> <li>-2 fixed compartments by 5 liters for developing and fixing</li> <li>-1 Fixed compartment by 28 lt for Intermediate and final washing</li> <li>- Automatic loading and overflow</li> </ul> <p><b>-W150-V220-Hz50</b></p> <p><b>DIMENSION :</b> Cm 27 x 48 x 96</p> <p><b>WEIGHT:</b> 12.5 Kg</p>	<p>Manual development of films rx up to the FORMAT CM 30 X 40</p> <p><b>INCLUDED:</b></p> <ul style="list-style-type: none"> <li>-2 Fixed compartments 2.5 lt for developing and fixing</li> <li>-1 Fixed compartment by 15 lt for intermediate and final washing</li> </ul> <p><b>DIMENSION :</b> Cm 22 x 38 x 53</p> <p><b>WEIGHT:</b> 11 Kg</p>	<p>Manual development of films rx up to the FORMAT CM 30 X 40</p> <p><b>INCLUDED:</b></p> <ul style="list-style-type: none"> <li>-2 Fixed compartments 2.5 lt for developing and fixing</li> <li>-1 Fixed compartment by 15 lt for intermediate and final washing.</li> </ul> <p><b>-W150-V220-Hz50</b></p> <p><b>DIMENSION :</b> Cm 22 x 38 x 53</p> <p><b>WEIGHT:</b> 11.5 Kg</p>

- Constituted PVC antacid, complete with removable lid for developing, fixing.
- Faucets antacid