

TECHNICAL DATA-SHEET

X-RAY FILM VIEWERS WITH SHUTTERS



- Stove enamelled steel structure, thickness 0,8 mm;
- Fire epoxy painting color peel RAL 9010;
- Smelted opaline plexiglass screen;
- Neon tubes of 6.000°Kelvin with high frequency reactors feeding;
- Upper and/or central film holder;
- Speed bipolar switches on the front complete with dimmer which allows a very good light regulation from 100% to 10%;
- External fuses lagging type 5x20 (positioned on the bottom right side);
- Feeding cable complete with schuko plug, CE 3. 0,75 ;
- Predisposed for wall installation through expansion dowels;
- Illuminator is made according to Directive 93/42/CEE – Medical Devices - Class 1
And Directive 73/23 CEE (Low tension) 89/336 CEE (Electromagnetic compatibility)
The construction meets Norms EN 60601 -1,2 / Din 6856 -1,2;
- Feeding 230V / 50Hz + earth;
- Maintenance : recurrent cleaning (each 6 months) of the internal screen.

X - RAY FILM VIEWERS WITH SHUTTERS

| CODICE | SIZE | OUTSIDE DIMENSION cm.. | WEIGHT Kg.. | NUMBER OF TUBES | LIGHT SECTION | W | A | FUSES | CD/MQ |
|---------------|-------------|---------------------------------------|------------------------|--------------------------------|--------------------------|----------|----------|--------------|--------------|
| LXJ40 | 40 X 43 | 60 X 62,5 X 10 | 20 | 4 – 15W | 1 – 40 X 43 | 65 | 0,35 | 1 | 4500 |
| LXJ80 | 80 X 43 | 100 X 62,5 X 10 | 31 | 4 – 30W | 1 – 80 X 43 | 135 | 0,60 | 1 | 6200 |
| LXJ120 | 105 X 43 | 125 X 62,5 X 10 | 38 | 4 – 36W | 1 – 105 X 43 | 165 | 0,75 | 1 | 5800 |
| LXJ121 | 120 X 43 | 140 X 62,5 X 10 | 43 | 4 – 36W | 1 – 120 X 43 | 165 | 0,75 | 1 | 5800 |
| LXJ140 | 140 X 43 | 160 X 62,5 X 10 | 51 | 4 – 58W | 1 – 140 X 43 | 225 | 1,10 | 2 | 7000 |

X – RAY FILM VIEWERS TO TWO SECTIONS COMBINED WITH OR WITHOUT SHUTTERS

Technical characteristics as described above, but consist of two superimposed negatoscopes, suitable to realize a single structure.

- Lower film viewer with 4 blackout curtains
- Superior film viewer without blackout curtains

| | | | | | | | | | |
|----------|----------|----------------|----|---------|-------------|-----|-----|---|------|
| LXJ140/2 | 144 X 83 | 168 X 134 X 12 | 70 | 8 – 58W | 2– 144 X 43 | 450 | 2.2 | 2 | 7000 |
|----------|----------|----------------|----|---------|-------------|-----|-----|---|------|