

TECHNICAL DATA-SHEET

X-RAY FILM VIEWERS WALL RECESSED SERIES



- Homogeneous lighting by n°3 fluorescent tubes with colour temperature 6000° K.
 - Pressing rolls film-holding system
 - Stove enamelled steel structure
 - Outside polished steel frame
 - White monochromatic diffuser
 - Equipped with terminal board for electrical connection.
 - Predisposed for explosion-proof installation.
 - Which this declaration refers, conforms to Directive 93/42/CEE – Medical Devices - Class 1
- Moreover meets the Directive 73/23 CEE (Low tension) 89/336 CEE (Electromagnetic Compatibility)

The construction conforms to Norms EN 60601 -1,2 - Din 6856 -1,2

TECHNICAL DATA - SHEET

WALL RECESSED SERIES

CODE	WALL HOLE NET LIGHT FRAME cm.	OUTSIDE FRAME DIM	NUMERO LAMPAD E	LIGHT SECTIONS	W	A	FUSES	CD/MQ
M-NI-90	Cm 95,5 x 50,5 x 12	Cm 102 X 57 X 11	3 – 30 W	1 – 90 X 43	120	0,6	1	3300
M-NI-120	Cm 125,5 x 50,5 x 12	Cm 132 X 57 X 11	3 – 36 W	1 – 120 X 43	130	0,6	1	3500