

CABLAS s.r.l.

Via Michelangelo Buonarroti, 52/B - 20093 Cologno M.se (Mi)

INTENSIFYNG SCREEN RARE EARTH BLUE / GREEN GREEN GRADUAL



INSTRUCTION AND INSTALLATION HANDBOOK



REV.01/2015

1.0 DESTINATION

The intensifying screens are intended to be inserted into x-ray cassettes for the radiographic analysis of human and animal body. Moreover are used in the industrial field.

2.0 CONSTRUCTION

The intensifying screens are compatible with all conventional x-ray films and are destined to improve the quality of overall rx image. Our range of screens is suitable for all type of x-ray cassettes.

Every screen provides high contrasts, details perception, even with low rx radiation dosage.

3.0 USE

Before using the x-ray cassette you must apply the intensifying screens inside the cassette (one side on the bottom of the cassette and one side under the cover).

Successively you must insert a virgin x-ray film.

It is recommended to avoid violent pushes that may deform and adulterate the normal use of the screens.

4.0 CLEANING

The screens do not need maintenance.

INTENSIFYING SCREEN RARE EARTH BLUE - SPEED 400

CODE	GND	N.REP.	SIZE
H-AS-01	Z11079008	304448	13 X 18
H-AS-02	Z11079008	304639	18 X 24
H-AS-03	Z11079008	304646	24 X 30
H-AS-04	Z11079008	304641	18 X 43
H-AS-05	Z11079008	304652	30 X 40
H-AS-06	Z11079008	304657	35 X 35
H-AS-07	Z11079008	304659	35 X 43
H-AS-08	Z11079008	304635	15 X 30
H-AS-09	Z11079008	304637	15 X 40
H-AS-10	Z11079008	304643	20 X 40
H-AS-12	Z11079008	304655	30 X 90
H-AS-13	Z11079008	304649	30 X 120

INTENSIFYING SCREEN RARE EARTH GREEN GRADUAL

CODE	CND	N REP.	SIZE
H-GR-09	Z11079008	304709	15 X 40
H-GR-10	Z11079008	304711	20 X 40
H-GR-12	Z11079008	304718	30 X 90
H-GR-13	Z11079008	304715	30 X 120
H-GR-14	Z11079008	304723	40 X 90
H-GR-15	Z11079008	304720	40 X 120

INTENSIFYING SCREEN RARE EARTH GREEN - SPEED 40

CODE	CND	N REP.	SIZE
H-AT-20	Z11079008	304661	13 X 18
H-AT-21	Z11079008	304674	18 X 24
H-AT-22	Z11079008	304681	35 X 43
H-AT-23	Z11079008	304676	18 X 43
H-AT-24	Z11079008	304683	30 X 40
H-AT-25	Z11079008	304698	35 X 35
H-AT-26	Z11079008	304704	35 X 43
H-AT-27	Z11079008	304669	15 X 30
H-AT-28	Z11079008	304671	15 X 40
H-AT-29	Z11079008	304672	20 X 40
H-AT-31	Z11079008	304693	30 X 90
H-AT-32	Z11079008	304684	30 X 120

CE CONFORMITY DECLARATION

The undersigned "CABLAS S.r.l."
office at Via Michelangelo Buonarroti 52/B
20093 Cologno Monzese (MI) - ITALY

DECLARES

Under own responsibility that the product

INTENSIFYING SCREENS RARE EARTH BLUE / GREEN - GRADUAL

Code

H-AT-20	H-AT-21	H-AT-22	H-AT-23	H-AT-24	H-AT-25
H-AT-26	H-AT-27	H-AT-28	H-AT-29	H-AT-31	H-AT-32
H-AS-01	H-AS-02	H-AS-03	H-AS-04	H-AS-05	H-AS-06
H-AS-07	H-AS-08	H-AS-09	H-AS-10	H-AS-12	H-AS-13
H-GR-09	H-GR-10	H-GR-12	H-GR-13	H-GR-14	H-GR-15

Which this declaration refers, conforms to
Directive 93/42/CEE – Medical Devices - Class 1

Moreover meets the Directive 73/23 CEE (Low tension)
89/336 CEE (Electromagnetic compatibility).

The construction conforms to Norms
EN 60601 -1,2 - Din 6856 -1,2

Cologno Monzese, year 2013

CABLAS S.R.L.
Via Michelangelo, 52/B
20093 COLOGNO MONZESE (MI)
S.r.l. - Via Sormani 34 - COLOGNO M.SE
Tel. P.A. 02 2538088 - Fax 02 27300644
Cod. Fisc. e P. IVA 11809230156



